

Steve Holland, numero uno CCD engineer at LBNL, will be giving the noon detector seminar at ANL, Thursday, April 23rd.


Will try to schedule a private meeting, next 2 slides are what I plan to show him...



Back-illuminated, left edge




Front-illuminated, left edge



Plan is to go back to APS with new front-illuminated CCD, and hopefully with a smaller beam

Received 6 CCDs from FNAL with unknown status, 5 are front-illuminated, first one tested looks funny




APS News


- Pinhole studies by Rich and Hal
- Continuation with <300um slit opening
- 2 Pinhole Plan?



APS News

- · Pinhole studies by Rich and Hal
- Continuation with <300um slit opening
- 2 Pinhole Plan?





- 3 options (?): 1) talk to Tim Graber about APS beamline 5um capability
 - 2) 1 pinhole fixed, 2nd on xy-stage (hard to reproduce)
 - 3) make tube for both pinholes (hard to point to beam spot)

Pinhole 1 radius	Pinhole 2 radius	Distance Pinhole 1 to CCD	Image radius on CCD	
10um	2.5um	15 cm	6um	
15um	2.5um	20 cm	6um	
30um	5um	30 cm	7um	

The Detroit Area Schools Supernova Search

26 Mar 2009

Adam Lincoln
David Cinabro
Matt Taylor
Nick Davis

Wayne State University
Dept. of Physics and Astronomy

DASSS: program overview

- Provide equipment to participating schools
 - 8 or 12 inch telescope, computerized mount, CCD
- Schools commit to weekly observation sessions for DASSS, but otherwise free to use equipment
- Once a week or so, a school takes images of about 10 assigned galaxies.
 - Each school pulls the target information from central website.

SK Comment: Based on discussion at meeting, this looked doomed to failure, will give up before seeing a SN. Seems much better to get to ~16 magnitude sensitivity and do SN follow-up based on IAU notification.



12" can get to 11.6 mag in 2 minutes, and 16" can get to 10.8 mag in 20 seconds (stars).


$$11.6 = Zeropoint12" + 2.5*log10(120s)$$
 Zeropoint12" = 6.4
Mag = 6.4 + 2.5*log10(3600s) = **15.3**

M51 taken with 16" and 3 filters



backup...

Copper Tube Spectrum








Figure 2.19 X-ray absorption length compared to visible and near-IR photons.

Front-illuminated



Back-illuminated, red=600nm, blue=900nm




ANL/HEP test-stand at APS x-ray lab.

Using x-ray tube plus 5 μm tungsten pinhole.


Energy range of x-rays gives similar absorption lengths as optical range (Janesick).


Sub-micron stepping motors move test-stand w/o disturbing beam.




CCD News

- Ordered 5 new pinholes (20,50,100,200,300 um) to help make sense of APS beamline, hope to put on motorized mounts
- Should get new front-illuminated CCDs from FNAL within a week, but need to test in HEP first.







JN 60

X-ray tube with standard focus focus size: 1 mm x 10 mm

X-ray tube with fine focus focus size: 0.4 mm x 8.0 mm

Anoc	de material	Tube output	RequisNo.	Anode material	Tube output	RequisNo.	
	Cu	2000 W	1900	Cu	1500 W	1907	
	W	2400 W	1901	W	2000 W	1908	
	Мо	2400 W	1902	Мо	2000 W	1909	
	Cr	1800 W	1903	Cr	1300 W	1910	
	Co .	1800 W	1904	Co c	1200 W	1911	
	Fe	1500 W	1905	Fe	900 W	1912	
	Ag	2000 W	1906	Ag	1500 W	1913	

Technical Data

Tube voltage:

maximum 60 kV

Anode angle:

zero degrees

Focus size at a radiation angle of 6°:

SN 60

Point focus

1.0 mm x 1.0 mm

Line focus

0.1 mm x 10.0 mm

SF 60

Point focus

0.4 mm x 0.8 mm

Line focus

0.04 mm x 8.0 mm

Inherent filtration:

0.4 mm Be

Filament current:

maximum 3.8 amps

(at approx. 10 V)

Weight

(SN 60 and SF 60): short anode 1.6 kilos

long anode 2.0 kilos

Water circulation


Cooling system: Cooling water

minimum 3.5 liters/min

requirement: Water intake


temperature:

maximum 35°C



Anode	Exciting voltage (kV)	Energy Κα ₁ (keV)	Kα, in Å	Κα ₂ in Å	Κā in Å	Kβ ₁ in Å	Beta Filter
Cr	6.00	5.41	2.28975	2.293660	2.2911	2.08492	V
Fe	7.11	6.40	1.936087	1.940025	1.9374	1.75665	Mn
Co	7.71	6.93	1.789007	1.792892	1.7903	1.62083	Fe
Cu	8.98	8.05	1.540598	1.544426	1.5418	1.392251	Ni
Мо	20.00	17.48	0.7093165	0.713607	0.7107	0.632303	Zr
Ag	25.52	22.16	0.5594205	0.563811	0.5609	0.497081	Rh
W _K	(69.51) (12.09)	59.31 8.40	1.4765	1.4875		1.2815	


Taken from: Int. Tab. of X-Ray Crystallography, Vol. IV, 1974 modified as with NBS, acc. to DESLATES, R.D., and HENIUS, A. (1973). Phys. Rev. Lett. 31, 972




CCD News


- APS analysis underway...
- Having 9 dead pixels is understood, 6 extra pixels in serial register, and Leach says other 3 are from the electronics pipeline. (latest CCD has 10 dead?)


4 pixel shift is global, due to extra 4 "prescan" pixels, easily flagged by glowing edge




Growing effective pixel size at edge, opposite of what we saw in front-illuminated








Beam Centroid From Encoder (Column)








CCD Quantum Efficiency Tests in HEP

- Front vs back-illuminated
- Hal's tests of integrating sphere throughput
- Hal's tests of monochromator passing "higher orders"

Juan says that front-illuminated QE is supposed to look like open triangles (which we measured), and very different than back-illuminated








backup

Precam dewar initial design

Vic Guarino (ANL), with input from Joe Bernstein, Herman Cease, Greg Derylo, Juan Estrada, Brook Gregory, ... has an initial 4-CCD dewar design.

Positronic connector for temperature readout and heating

Slots based on DES V2 focal plane support plate



CCD TODO Rest of 2009...

- Writeup Edge Effects in journal article (need more data? Depends on LBL reaction to draft)
- Meet with APS group on future work, new ~50-100 um pinhole on moving stand instead of slits (+ 5um pinhole)? Goal easy to achieve <=10um beam.
- Test new pinhole setup with Weizeorick FastCCD or APS PMT (rates only) or both
- Get new front-illuminated CCD from FNAL
- HEP teststand: QE mystery, DES overflow, modify focal plane for science packages, rest of FNAL testing procedure
- Precam

Precam dewar initial design

Plan:

- 1) Discuss more with experts (for example, 0.5" clearance for CCD cables to make 90 degree bend?)
- Make the focal plate in March and test it in current ANL dewars
- 3) Build the dewar this summer independent of precam decision (already funded for FY09)
- 4) Take to CTIO in December as part I of phased approach? (if Precam approved...)