

# System Optimization for Time-of-Flight PET

**William W. Moses**

*Lawrence Berkeley National Laboratory*

February 28, 2009


## Outline:

---

- PET Electronics Requirements
- PET Electronics Trends
- OpenPET Electronics

- This work was supported in part by the U.S. DOE (contract No. DE-AC02-05CH11231) and in part by the NIH (NIBIB grant No. R01-EB006085).
- Thanks to S. Choong, M. Ullisch, M. Janecek, and C. Vu of LBNL, C. Melcher of U. Tennessee, P. Szupryczynski and M. Andreaco of Siemens Medical Solutions, and S. Hales of Hamamatsu Photonics.

# Time-of-Flight in HEP




- Photons Produced Promptly
- Photons Travel in ~Same Direction
- Small Time Variations due to Path Length Difference
- Small Variations due to Photon Production Position

**Time Spread Between Photons  
Arriving at PMT is Small**


# Arrival Time at PMT in HEP

(4 Photons Produced)



- **Within Event, Photons Usually Overlap in Time**
- **Small Event-to-Event Variations**

# Time-of-Flight in PET




- Photons Produced w/ Scintillator Decay Time
- Photons Travel in All Directions
- Large Time Variations due to Path Length Difference
- Large Variations due to Photon Production Position

**Time Spread Between Photons  
Arriving at PMT is Large**

# Arrival Time at PMT in PET


(4 Photons Produced)



- **Within Event, Photons Often Don't Overlap in Time**
  - **Large Event-to-Event Variations**

# Why Good Timing Is Hard in PET

Tek **Stop:** 2.00GS/s



**Statistical Fluctuations Make Accurate Timing Difficult!**

# Measure Each Component in Timing Chain



- Construct a Chain with “High Performance” Components
- Measure “High Performance” Timing Resolution
- Replace One of the “High Performance” Components with a “Production” Component
- Re-Measure Timing Resolution
- Difference (in Quadrature) is the Timing Contribution of the “Production” Component

**Identify Where Can We Make Improvements**

# “Reference” Timing Signal

Ortec 556 TAC  
16-bit ADC



Tennelec /  
Canberra 454




Hamamatsu H-5321


140 ps TTJ  
700 ps Rise Time  
1" Quartz Window



1x1x1 cm BaF<sub>2</sub>



<sup>68</sup>Ge





# “High Performance” Measurement

3x8x12 mm  
LSO


Hamamatsu H-5321

140 ps TTJ  
700 ps Rise Time


1" Quartz Window



Tennelec /  
Canberra 454



Ortec 556 TAC  
16-bit ADC




# Scintillator Crystal

Hamamatsu H-5321

Tennelec /  
Canberra 454

Ortec 556 TAC  
16-bit ADC

6.75 x 6.75 x 25 mm LSO



394 ps – 221 ps  $\Rightarrow$   
**326 ps fwhm (batch 1)**

502 ps – 221 ps  $\Rightarrow$   
**451 ps fwhm (batch 2)**

# Scintillator “Block” Effects

**CPS ACCEL V1 Block**

**8x8 Array**


**6.75 x 6.75 x 25 mm<sup>3</sup> LSO**

**Hamamatsu Fast 2” PMT**

**250 ps TTJ**

**1500 ps Rise Time**

**Prototype On Loan**



**Configuration**

**Single Crystal  
Block**

**Resolution**

**355 ps**

**577 ps**

**⇒ Block Contribution**


**454 ps**

# What Limits Timing Resolution?

## Non-TOF Block Detector Module

Baseline  
160 ps

Crystal  
Geometry  
326 ps




Light Sharing  
454 ps

PMT  
422 ps

PMT Array  
274 ps


- Many Factors
- “Optical Geometry” Particularly Important

# Proposed Side-Coupled Design



**Shorter Optical Path Length & Fewer Reflections**


# Optimization: LSO Composition



- Both Scintillators Have Same Light Output (photons/MeV)
- Red Decay Time is 2x Longer Than Blue Decay Time

- Predicted Timing Resolution  $\propto 1/\sqrt{I_0}$
- Want High Total Light Output & Short Decay Time
- Possible By Co-Doping LSO With Calcium

# Optimization: LSO Composition



- **Ca-Doping Gives Good Timing Resolution**
  - **~15% Improvement Over Normal LSO**




# Optimization: Photomultiplier Tube

## Blue Sensitivity Index

## R9800 Spectral Response


Peak QE



- Predicted Timing Resolution  $\propto 1/\sqrt{\text{QE}}$
- Want High Quantum Efficiency Version of PMT



# Measured Results: High QE PMTs



- Increased QE Improves Timing Resolution by 7%
- Expect 10% Improvement with 35% SBA PMT

# Summary

Hardware	Single (ps fwhm)	Coinc. (ps fwhm)	TOF Gain
End-Coupled Crystal	384	544	4.3
Side-Coupled Crystal	218	309	7.6
Co-Doped LSO	182	258	9.1
32% QE PMT	155	219	10.6
35% QE “SBA” PMT	148	209	11.1

- TOF PET with *Significantly* Better Timing is Possible
  - To Achieve, We Must Optimize the *System*

# TOF Gain for Whole-Body PET (35 cm)

Hardware	$\Delta t$ (ps)	TOF Gain
BGO Block Detector	3000	0.8
LSO Block (non-TOF)	1400	1.7
LSO Block (TOF)	550	4.2
LaBr <sub>3</sub> Block	350	6.7
LSO Single Crystal	210	11.1
LuI <sub>3</sub> Single Crystal	125	18.7
LaBr <sub>3</sub> Single Crystal	70	33.3


***Incredible Gains Predicted...***

# Detector Module Design

**Two LSO Crystals**  
(each  $6.15 \times 6.15 \times 25 \text{ mm}^3$ )

**Reflector**  
(on all five faces of  
each crystal, including  
the face between the  
two crystals)

**Optical Glue**  
(between lower crystal  
faces and PMT)



**Hole in Reflector**  
**On Top Face of**  
**Crystals**


**PMT**  
(Hamamatsu  
R-9800)

**Two Side-Coupled Scintillator Crystals per PMT**

# Detector Ring Geometry

**Crystal of Interaction**


**Exploded View**



- Top face of each crystal (with hole in reflector) is coupled via a small ( $<1$  mm) air gap to the edge of one opposing PMT.
- Light seen by the opposing PMT is used to decode the crystal of interaction.

**Crystals Decoded by Opposing PMT**


# Camera Geometry



- Detector ring is 825 mm diameter, 6.15 mm axial
- 192 detector modules, 384 LSO scintillator crystals
- Adjustable gap (6 – 150 mm) between lead shields allows “scatter-free” and “3-D” shielding geometries


**“Real” Single-Ring PET Camera for Humans & Phantoms**

# Timing Values will be for a *Single* Detector\*



- Only Accept Events in Photopeak Window
- Subtract (in Quadrature) 150 ps Trigger Contribution

# Light Transport Affects Timing Resolution



**Long, Thin Crystals Have Slower Rise Time**