

# Towards an integrated high resolution TOF PET system

Jean-Francois Genat

University of Chicago


February 28<sup>th</sup> 2009

# Pico-second timing and TOF PET

Pico-second timing of FWHM  $\Delta t$  can provide:

Coordinate along the 2-photons : FWHM improved by  $2D/c\Delta t$

Transverse position  
using T Lines



Intrinsic resolution (FWHM)  
e<sup>+</sup> / e<sup>-</sup> 1mm ( $\Delta t=66\text{ps}$ )  
 $\delta \Theta$  4mR

- Improve image quality,
- Speed up the convergence of the reconstructing algorithms

Best resolution of **500-300ps** FWHM, 7.5-4.5cm with LSO or LaBr<sub>3</sub> crystals  
<1mm transverse (T-lines), centroids in the axial direction

# TOF vs No-TOF PET

- Surti et al. NSS MIC 2006




Fig. 6. Transverse slices for the reconstructed images of the 35-cm diameter lesion phantom with Non TOF (A) and TOF (B) reconstructions.

$$\text{TOF: } \Delta x = c \Delta t / 2$$

$\Delta t$  due mainly to **crystal**

# Recent progress

## **INNOTEP Project**

**G Montarou**

### **Monte Carlo assessment of time-of-flight benefits on the LYSO-based discovery with PET/CT scanner**

Geramifar, P. et al

Biomedical Imaging: From Nano to Macro, 2008. Volume , Issue , 14-17 May 2008 pp364-367

### **Potential Advantages of Digitally Sampling Scintillation Pulses in Timing Determination in PET**


Xie Quingguo et al

IEEE Trans on Nucl Sci 2008


# Crystals signals

Block/Single Crystal @ 511 KeV (%)	Timing res. (ps)	Light output $10^4$ ph/MeV	Energy res.
---------------------------------------	------------------	-------------------------------	-------------

- |                     |        |     |    |
|---------------------|--------|-----|----|
| • BGO               | 3000   | 0.8 | 12 |
| • LSO               | 550    | 2.5 | 10 |
| • LuI <sub>3</sub>  | 125    | 9   | 6  |
| • LaBr <sub>3</sub> | 350/70 | 6   | 3  |
- W. Moses, Oct 2008



LSO+PMT rise time is ~1.4ns



LSO only is 1ns

H Kim,  
Chicago

# Timing and position resolution limits for PET

## -1 Along the $\gamma$ ray axis

Limited by crystals timing resolution

$\text{LuI}_3$  : 125ps resolution with 3 times more photons than LSO, but are too expensive.

$\text{LaBr}_3$  provides 350ps resolution at 2/3 less photons, good trade-off, compared to LSO

350ps is **5.2cm FWHM**

## -2 Transverse

Micro-Channel Plates + transmission lines read with sampling chips in the 10GHz range can provide 10ps timing (laser), so transverse position measurements to  $500\mu\text{ m}$  FWHM, assuming T-lines velocity  $c/3$ , and crystal resolution not contributing since two ends signals are correlated. Crystal size quantization to be added:

total is  $650 + 500 = 820\mu\text{ m}$  FWHM

## -3 Axial

Crystal cannot be smaller than  $1 \times 1\text{mm}^2$  ( $650\mu\text{ m}$  FWHM)

# Fast Photo-detectors


30PEs	Time (FWHM)	Position (FWHM)	Area	Pix size mm <sup>2</sup>
(Ma)PMT	60ps	9.1mm	2" x 2"	6 x 6
SiPM	70ps	10.6mm	linear array, 32ch	1 x 1
MCP	25ps	3.8mm	2" x 2"	1.6 x 1.6

MCP best for large size, timing and pixel size

# Position + time with Transmission Lines

- Pico-second timing provides
  - **Fast timing**
  - **Transmission lines readout.**
    - Transverse position along the line**
    - Axial coordinate from centroids**

Less electronics channels for large area sensors



$$\frac{1}{2} (t_1 + t_2) = \text{time}$$


$$t_1 - t_2 / t_1 + t_2 = \text{transverse position}$$

$$\Sigma \alpha_i a_i / \Sigma \alpha_i = \text{axial position}$$




# MCP + TLs

Laser @ 50PEs:    70ps FWHM one-side    10ps diff




# Timing with laser on 2" x 2" MCPs (TLs on Printed Circuit Board)

Fitting Time resolution v. Light for MCP with Delay Line Readout





# T-lines two-ends with 25 $\mu$ m MCP




25  $\mu$  m pore MCP signal at the output of a ceramic transmission line  
Laser 408nm, 50 $\Omega$  , no amplification

# Signal processing for pico-second timing


- Multi-threshold
- Constant Fraction
- Sampling and Waveform analysis




# Pulse Sampling and Waveform Analysis



# MCP + TL (laser)




25um tube, 50 PEs  
FWHM = 10ps      500  $\mu$  m




10um tube  
FWHM = 8.9ps      450  $\mu$  m

**NB. Not obtained from crystal signals, but laser fast light source**

# A MCP based Micro-PET system (H. Frisch)



# A PET system with on-detector electronics





# Self-triggered Readout Architecture

Strategy: Extract from the front-end relevant data only

- Sample and overwrite continuously in the 10 GS/s range (analog)
- Digitize in the front-end any activity above a given threshold
  - Using front-end sampling chip triggering capability  
and *optionally*
  - Process on front-end with FPGAs to extract time, amplitude
- Send to system crate for processing via a few optical/*copper* links
- Process data in ATCA System crate
  - Coincidences, reconstruction, calibrations
- Timing distribution from system crate needed with a  $\sim 10$ ps stability
  - Not too critical on a relatively small system

# Sampling Chips

State of the art:

	Sampling	Bandwidth	Dyn. range	Depth	PLL	ADC	Trigger
	GHz	GHz	bits			bits	
G. Varner (Hawaii)	6	1.0	10	1024	no	12	experience
S. Ritt (PSI)	6	.8	11.5	256	3.9ps		no
D. Breton (Orsay) + E. Delagnes (Saclay)	2.5	.3	13.4	250	20ps	no	no

# 130nm CMOS Timing Sampling chip

- 40 GHz sampling rate
- Analog bandwidth ~ 1.5 GHz
- Triggered mode ←
- Depth 256
- On-chip 10-bit ADC ←

Presently under development at U Chicago

# Dataflow

## Dataflow

- Assume 64 samples (taken for 1.6 ns at 20GS/s), in nine adjacent crystals
- Assume 1MCps rate, a safety factor of four (noise + system noise)
- Each of four sections linked to one  $\mu$  TCA crate runs at:



Without sample processing

$$4\text{M Cps} \times 2 \times 4 \text{ layers} \times 9 \text{ crystals} \times 64 \text{ samples} \times 8\text{bits} = 72 \text{ Gb/s}$$

With sample processing (FPGA)

$$4\text{M Cps} \times 4 \text{ layers} \times 9 \text{ crystals} \times 4 \text{ words} \times 8\text{bits} = 4.6 \text{ Gb/s}$$

Deadtime



# ATCA System

- **Use ATCA Advanced Telecom Computing Architecture**
  - ATCA with the full sampled data readout option ( $>70$  Gb/s input)
  - Use fibers at 10 Gb/s (Gigabit Ethernet)
- Micro-TCA system may be sufficient in the case of on-detector digital processing ( $> 4$  G/s input)
- A few serial frast cables sufficient

High resolution image output

# Conclusion

Use of an integrated PET system making use of Micro-channel Plates sensors, transmission lines readout, custom sampling electronics, and digital signal processing allows :

- Extracting timing to sub-millimeter precision in the transverse dimension with a reduced number of electronics channel,
- Reducing the amount of data to be read and processed off-line by one order of magnitude at least.

# Extra slides

# 10um tube

- Tuning the HV at 2.5 kV to have the same MCP pulse amplitude as the 25um tube at 2kV (30mV)

50PEs, 2.5kV  
std=3.82ps  
vs 2.5ps (simulation)

18PEs, 2.5 kV  
std = 6.05ps  
vs 7ps (simulation)

